

# MEDIA RELEASE



|                          |  |            |         |
|--------------------------|--|------------|---------|
| <b>Release Status</b>    | For immediate release  | <b>Ref</b> | 1502DIA |
| <b>Product Enquiries</b> | <b>Diatron Group</b><br><b>International contact:</b><br>Brian Conibere<br>International Sales and Marketing Manager<br>t: +44 1234 880697                      m: +44 7717 856815<br>e: <a href="mailto:brian.conibere@diatron.com">brian.conibere@diatron.com</a>  |            |         |
| <b>Media enquiries</b>   | <b>Alto Marketing Limited</b><br>Victoria Coupe<br>t: +44 (0)1489 557672                      e: <a href="mailto:victoriac@alto-marketing.com">victoriac@alto-marketing.com</a><br>Twitter: <a href="https://twitter.com/AltoMarketing">@AltoMarketing</a> w: <a href="http://www.alto-marketing.com">www.alto-marketing.com</a> |            |         |

## Diatron to launch new 5-part veterinary hematology analyzer at Medica

The Abacus Vet5 provides excellent usability and LIS integration



[Contact [victoriac@alto-marketing.com](mailto:victoriac@alto-marketing.com) for high res image]

**BUDAPEST, HUNGARY (November 2, 2015):** The [Diatron Group](http://www.diatron.com), global provider of hematology and clinical chemistry analyzers, will launch a new 5-part differential veterinary hematology analyzer - the [Abacus Vet5](http://www.diatron.com), in Dusseldorf at the forthcoming Medica show on the 16-19 November. This will be introduced on **booth 3A 07** alongside Diatron's new 3-part differential human hematology analyzer, the [Diatron Aquila](http://www.diatron.com), which delivers lower costs and wider range of use, from small laboratory to point-of-care environments.

The Abacus Vet5 analyzer has a completely new look and is operated via a large capacitive touch screen, greatly enhancing user experience. New features include: increased on-board data storage capability of up to 5000 records including histograms;

# MEDIA RELEASE



a new comprehensive QC package; and more data export/connectivity functionalities, providing improved communication with laboratory information systems (LIS).

The new analyzer uses well-proven Diatron Veterinary reagents which enable a wide range of hematological parameters to be tested on a broad panel of animal species. Each reagent pack for the new Abacus Vet5 includes a special hardware key which transfers important product data (e.g. lot number and expiry date) automatically to the analyzer, eliminating transcription errors while protecting the system from being used with non-recommended reagents.

### *New 3-part differential human hematology analyzer - Diatron Aquila*

Also on booth 3A 07, the **Diatron Aquila**, a [new 3-part differential hematology analyzer](#) (60 test /hour throughput) will be introduced to the laboratory world for the first time. This innovative analyzer has been specially designed to provide unique solutions to the many needs of the rapidly evolving worldwide hematology markets. Requiring minimal bench space, this small analyzer has low reagent consumption (<10ml), needs only a small sample size (10ul) and has an on-board reagent pack\*\* (sufficient for circa 200 measurements).

For more information, please visit booth 3A 07 at Medica or the Diatron Group website: [www.diatron.com](http://www.diatron.com).

\*\* *patent applied for on the on-board reagent pack*

## **Editors' Notes**

### **The Diatron Group**

Diatron develops, manufactures and markets hematology and clinical chemistry analyzers, and associated reagents for human medical and veterinary use. The company was founded over 20 years ago and, since its inception, has been at the forefront of laboratory diagnostics. Today, in its new production facility in Budapest, Hungary, the extensive R&D/Engineering team develops and manufactures a broad product range including: hematology analyzers, hematology reagents (both for its own and other manufacturers' analyzers), hematology control material, clinical chemistry analyzers, clinical chemistry reagents and clinical chemistry controls.

The brand name of Diatron has been established throughout the world with its products being sold and marketed in more than 100 countries. Today, there are over 30,000 Diatron clinical chemistry and hematology analyzers in laboratory use, with a customer base that continues to grow strongly year after year.

# MEDIA RELEASE



Address: Tablas u. 39. , H-1097 Budapest, Hungary

Phone: +36 1 436 9800

Email: [sales@diatron.com](mailto:sales@diatron.com)

DIATRON USA

Address: 4548 N. Hiatus Rd., Sunrise, FL 33351, USA

Phone: 877-684-1139

Email: [sales@diatron.com](mailto:sales@diatron.com)