

Diatron presents the first in a new family of hematology analyzers – the Aquarius 3 Three-part differential

Brings unrivalled simplicity plus renowned efficiency to the decentralized market segment.

Simply counts for you!



Diatron announces the launch of the Aquarius 3, its latest IVD analyzer for Complete Blood Count testing. In line with its predecessors, the new Aquarius 3 is easy to use, flexible and suitable for a variety of uses such as in doctor's offices, low to mid-size laboratories and other decentralized healthcare facilities.

Evolving from the Abacus 380, Diatron's flagship analyzer, the Aquarius 3 incorporates similar robust and trusted counting technology but is equipped with new features to improve operability and connectivity allowing it to adapt to diverse environments for the benefit of both customers and users alike.

While still retaining the important features of the Abacus 380, the Aquarius 3 has many other product enhancements, with the most obvious change being the addition of an 8" high resolution, color touch screen, making it clear and simple to operate. Others include

- 1D/2D barcode scanner for sample ID, QC and reagent management (closed system),

AQUARIUS 3 PRESS RELEASE

- 25 µl sampling volume, including micro-tube capability,
- Increased database capacity: 10,000 + USB mass storage,
- Wider choice of external printers

About Diatron MI

Located in Budapest, Hungary, Diatron is a member of the Stratec SE Group, Birkenfeld, Germany. It offers the complete package ranging from innovative R & D, excellence in manufacturing and reagent production plus the highest quality sales office and service support in the hematology, clinical chemistry and CLIA fields. This is reflected with products present in more than 110 countries distributed through an established and dedicated distributor network. With over 50,000 analyzers in laboratory use today, Diatron has a large customer base which is continuously growing year on year.

Contact

DIATRON INTERNATIONAL OFFICE

Address: Tablas u. 39, H-1097 Budapest, Hungary

Phone: +36 1 279 4400

Email: sales@diatron.com

Website: www.Diatron.com